

## Product Description

Product Name: LL/2-Fluc-Neo/eGFP-Puro  
 Catalog Number: CL073  
 Lot Number: CL-IM68

Species: Mouse (*Mus musculus*)  
 Strain: C57BL  
 Cell type: Lewis lung carcinoma  
 Parental cells: LL/2 (ATCC® CRL-1642™)  
 Morphology: Epithelial  
 Growth mode: Loosely adherent  
 Reporter genes: Firefly luciferase (Fluc)  
 Enhanced green fluorescent protein (eGFP)  
 Selection genes: Neomycin (Neo)  
 Puromycin (Puro)

This is a polyclonal population derived from the Lewis lung carcinoma LL/2 cell line (ATCC® CRL-1642™). Parental LL/2 cells were transduced with 1) LV-Fluc-P2A-Neo (Imanis #LV011) encoding the firefly luciferase (Fluc) cDNA under the spleen focus-forming virus (SFFV) promoter linked to the neomycin resistance gene (Neo) by a P2A cleavage peptide and 2) LV-eGFP-PGK-Puro (Imanis #LV031) encoding the enhanced green fluorescent protein (eGFP) cDNA under the SFFV promoter and the puromycin resistance gene (Puro) under the phosphoglycerate kinase (PGK) promoter. High Fluc and eGFP expressing cells were selected using G418 and puromycin. The lentiviral vectors are self-inactivating (SIN) vectors in which the viral enhancer and promoter have been deleted. Transcription inactivation of the LTR in the SIN provirus increases biosafety by preventing mobilization by replication competent viruses and enables regulated expression of the genes from the internal promoters without *cis*-acting effects of the LTR<sup>1</sup>.

## Mycoplasma Testing

The LL/2-Fluc-Neo/eGFP-Puro cell line has been tested for mycoplasma contamination and is certified mycoplasma free.

## Cell Line Authentication

In light of studies suggesting that 18-36% of cell lines utilized in biomedical research are contaminated or completely misidentified,<sup>2,3</sup> several funding organizations, including NIH, as well as major publishers, including those affiliated with the American Associate for Cancer Research (AACR), require cell lines used in research to be authenticated prior to publication<sup>4,5</sup>. The parental LL/2 cell line used to generate LL/2-Fluc-Neo/eGFP-Puro was purchased directly from ATCC.

## Recommended Uses

LL/2-Fluc-Neo/eGFP-Puro cells are suitable for *in vitro* and *in vivo* experimentation.

LL/2 cells form tumors post implantation into syngenic C57BL mice<sup>6</sup> and can be used for studying metastasis and cancer chemotherapeutics. The Fluc transgene in the LL/2-Fluc-Neo/eGFP-Puro cells facilitates noninvasive bioluminescent imaging of implanted cells, while the eGFP gene facilitates conventional microscopy of excised tissue.

## References

- <sup>1</sup>Miyoshi et al. J Virol. 1998. 72:8150-8157.
- <sup>2</sup>Hughes et al. BioTechniques 2007. 43: 575-586.
- <sup>3</sup>Chatterjee et al. Science 2007. 315:928-931.
- <sup>4</sup><https://grants.nih.gov/grants/guide/notice-files/NOT-OD-08-017.html>
- <sup>5</sup><http://www.aacrjournals.org/site/InstrAuthors/fora.xhtml#celllineuse>
- <sup>6</sup>Bertram et al. Cancer Letters, 1980. 11:63-73.

## Storage Instructions

Remove cells from the dry ice packaging and immediately store in the vapor phase above liquid nitrogen (below -130°C).

## Complete Growth Medium

Dulbecco's Modified Eagle's Medium (DMEM)  
 10% fetal bovine serum (FBS)  
 1% Penicillin/Streptomycin  
 1.25 mg/mL G418 (to maintain high Fluc expression)  
 2 µg/mL puromycin (to maintain high eGFP expression)

## Thawing Instructions

1. Thaw cells by gently swirling in a 37°C water bath. To limit contamination, do not submerge the O-ring and cap.
2. When cells are ~70% thawed (less than 1 min), remove the vial and wipe down with 70% ethanol. Allow tube to dry completely.
3. In a biosafety cabinet, transfer the cells into a 15 mL conical tube containing 5 mL of pre-warmed complete growth medium without selection drugs. Centrifuge at ~250 x g for 3-5 min.
4. Remove supernatant and resuspend cells in 1 mL complete growth medium without selection drugs. Transfer cells to a T75 flask containing 10 mL pre-warmed complete growth medium without selection drugs.
5. Incubate the culture at 37°C with 5% CO<sub>2</sub>. After 48 hours, replace the culture supernatant with complete growth medium containing 2 µg/mL puromycin and 1.25 mg/mL G418. Cells should reach full confluency 3-4 days after thawing.

## Subculturing Instructions

Volumes are given for a T75 flask. Increase or decrease as needed. In order to maintain high Fluc and eGFP expression, it is recommended that cells be subcultured in the presence of 2 µg/mL puromycin and 1.25 mg/mL G418. LL/2-Fluc-Neo/eGFP-Puro cells are loosely adherent, and may begin to detach from flasks prior to reaching 100% confluency. Therefore, cells should be passaged when they reach ~80% confluency.

1. Remove culture medium from cells.
2. Carefully wash the cell monolayer with 5-10 mL of phosphate buffered saline.
3. Add 2 mL of 0.25% Trypsin-EDTA solution to the flask and incubate at 37°C until cells have dissociated (approx. 2-5 min).
4. Neutralize the trypsin by adding 8 mL complete growth medium, and mix by gently pipetting up and down.
5. Transfer desired portion of the cells to a fresh T75 flask. Add fresh complete growth medium to a total volume of 10 mL and return cells to 37°C/5% CO<sub>2</sub> incubator.

For maintenance a subcultivation ratio of 1:10 is recommended. At this ratio cells will be ready for passage every 3-4 days.

## Freezing Medium

LL/2-Fluc-Neo/eGFP-Puro cells can be amplified and used to generate additional frozen stocks. Frozen stocks should be preserved in a designated cryopreservation medium or in complete growth medium without selection drugs supplemented with 5-10% DMSO.

## Additional Considerations

LL/2-Fluc-Neo/eGFP-Puro cells are loosely adherent. Coating culture plates with poly-D-Lysine prior to use can be used to increase cell adherence to the plates if necessary.

## Certificate of Analysis


Testing performed by Imanis Life Sciences

Test description	Result
Post thaw viable cell recovery	Pass QC
Sterility	No contamination detected
Mycoplasma	No contamination detected
Neomycin selection	Pass QC
Puromycin selection	Pass QC
Luciferase expression	Pass QC
Fluorescence expression	Pass QC

## Morphology:


Low density, 20X

High density, 20X




Low and high density photos taken 24 and 72 hours after thawing, respectively.

## Luciferase Expression



$10^4$ ,  $10^5$ , or  $10^6$  cells were placed in wells of a 96-well plate and 0.3 mg of d-luciferin was added to the indicated wells. The plate was immediately imaged using a Xenogen IVIS Spectrum.

## Fluorescence Expression



LL/2-mNIS-Neo/eGFP-Puro (green) or isotype control (LL/2-Fluc-Puro; grey) cells were fixed with paraformaldehyde and analyzed by flow cytometry (20,000 events).

## Legal Disclaimers

### Limited Product Warranty

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Imanis. Imanis shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

For *in vitro* use only. This certificate is a declaration of analysis at the time of manufacture.

### PURCHASER NOTIFICATION

#### LIMITED LICENSE NOTICE – RESEARCH USE ONLY

IMANIS LIFE SCIENCES HAS A LIMITED LICENSE UNDER PATENTS OWNED BY THE SALK INSTITUTE FOR BIOLOGICAL STUDIES THAT PERMITS IMANIS LIFE SCIENCES TO SELL PRODUCTS CONTAINING WPRE FOR RESEARCH USE ONLY AND NOT FOR ANY COMMERCIAL USES. EXCLUDED COMMERCIAL USES INCLUDE WITHOUT LIMITATION MANUFACTURING, PROVIDING A SERVICE, THERAPEUTIC, DIAGNOSTIC AND PROPHYLACTIC USES, AND ANY OTHER COMMERCIAL USES. USE OF THIS PRODUCT BY A PURCHASER FOR ANY PURPOSE OTHER THAN FOR RESEARCH IS UNAUTHORIZED AND PROHIBITED.

The Salk Institute actively licenses its patents for commercial uses, and a commercial use license may be available for Salk's WPRE patents. If you wish to inquire about such a license, please contact:

Office of Technology Development  
The Salk Institute for Biological Studies  
10010 North Torrey Pines Road  
La Jolla, CA 92037  
Phone: (858) 453-4100 extension 1278  
Fax: (858) 546-8093

Quality control by: RLV  
Quality Assurance by: SPR  
Effective Date: 11/25/15