# CT26.WT-iRFP-Neo



# **Product Description**

Product Name: CT26.WT-iRFP-Neo

Catalog Number: CL091 Lot Number: CL-IM88

Species: Mouse (Mus musculus)

Strain: BALB/c

Cell type: Colorectal carcinoma

Parental cells: CT26.WT (ATCC® CRL-2638™)

Morphology: Epithelial Growth mode: Adherent

Reporter gene: Near infrared fluorescent protein (iRFP)

Selection gene: Neomycin (Neo)

This is a polyclonal population derived from the murine colorectal carcinoma CT26.WT cell line (ATCC® CRL-2638<sup>TM</sup>). Parental CT26.WT cells were transduced with LV-iRFP-P2A-Neo (Imanis #LV033) encoding the near infrared fluorescent protein (iRFP; ex/em = 690/710 nm) cDNA under the spleen focus-forming virus (SFFV) promoter linked to the neomycin resistance gene (Neo) via a P2A cleavage peptide. High iRFP expressing cells were selected using G418. The lentiviral vectors are self-inactivating (SIN) vectors in which the viral enhancer and promoter have been deleted. Transcription inactivation of the LTR in the SIN provirus increases biosafety by preventing mobilization by replication competent viruses and enables regulated expression of the genes from the internal promoters without *cis*-acting effects of the LTR<sup>1</sup>.

# **Mycoplasma Testing**

The CT26.WT-iRFP-Neo cell line has been tested for mycoplasma contamination and is certified mycoplasma free.

### **Cell Line Authentication**

In light of studies suggesting that 18-36% of cell lines utilized in biomedical research are contaminated or completely misidentified, 2.3 several funding organizations, including NIH, as well as major publishers, including those affiliated with the American Associate for Cancer Research (AACR), require cell lines used in research to be authenticated prior to publication 4.5. The parental CT26.WT cell line used to generate CT26.WT-iRFP-Neo was authenticated and certified free of interspecies cross contamination by STR profiling with 27 STR loci.

### Recommended Uses

CT26.WT-iRFP-Neo cells are suitable for *in vitro* and *in vivo* experimentation.

CT26.WT cells form tumors and pulmonary metastases post implantation into syngenic BALB/c mice<sup>6</sup>. The iRFP transgene in the CT26.WT-iRFP-Neo cells facilitates *in vivo* noninvasive imaging of implanted cells.

# References

<sup>1</sup>Miyoshi et al. J Virol. 1998. 72:8150-8157.

<sup>2</sup>Hughes et al. BioTechniques 2007. 43: 575-586.

<sup>3</sup>Chatterjee et al. Science 2007. 315:928-931.

4https://grants.nih.gov/grants/guide/notice-files/NOT-OD-08-017.html

5http://www.aacrjournals.org/site/InstrAuthors/ifora.xhtml#celllineuse

<sup>6</sup>Wang et al. J Immunol, 1995. 154:4685-4692.

## **Storage Instructions**

Remove cells from the dry ice packaging and immediately store in the vapor phase above liquid nitrogen (below -130°C).

## **Complete Growth Medium**

Dulbecco's Modified Eagle's Medium (DMEM) 10% fetal bovine serum (FBS) 1% Penicillin/Streptomycin 0.4 mg/mL G418 (to maintain high iRFP expression)

# **Thawing Instructions**

- 1. Thaw cells by gently swirling in a 37°C water bath. To limit contamination, do not submerge the O-ring and cap.
- 2. When cells are ~70% thawed (less than 1 min), remove the vial and wipe down with 70% ethanol. Allow tube to dry completely.
- 3. In a biosafety cabinet, transfer the cells into a 15 mL conical tube containing 5 mL of pre-warmed complete growth medium without G418. Centrifuge cells at ~250 x g for 3-5 min.
- Remove supernatant and resuspend cells in 1 mL complete growth medium <u>without G418</u>. Transfer cells to a T75 flask containing 10 mL pre-warmed complete growth medium without G418.
- Incubate the culture at 37°C with 5% CO<sub>2</sub>. After 48 hours, replace the culture supernatant with complete growth medium containing 0.4 mg/mL G418. Cells should reach full confluency 2-3 days after thawing.

# **Subculturing Instructions**

Volumes are given for a T75 flask. Increase or decrease as needed. In order to maintain high iRFP expression, it is recommended that cells be subcultured in the presence of 0.4 mg/mL G418. CT26.WT-iRFP-Neo cells do not form a 100% confluent monolayer. It is recommended that cells be passaged when they reach ~80% confluency.

- 1. Remove culture medium from cells.
- Carefully wash the cell monolayer with 5-10 mL of phosphate buffered saline.
- 3. Add 2 mL of 0.25% Trypsin-EDTA solution to the flask and incubate at 37°C until cells have dissociated (approx. 1-2 min).
- 4. Neutralize the trypsin by adding 8 mL complete growth medium, and mix by gently pipetting up and down.
- 5. Transfer desired portion of the cells to a fresh T75 flask. Add fresh complete growth medium to a total volume of 10 mL and return cells to 37°C/5% CO<sub>2</sub> incubator.

For maintenance a subcultivation ratio of 1:10 is recommended. At this ratio cells will be ready for passage every 3-4 days.

# **Freezing Medium**

CT26.WT-iRFP-Neo cells can be amplified and used to generate additional frozen stocks. Frozen stocks should be preserved in a designated cryopreservation medium or in complete growth medium without G418 supplemented with 5-10% DMSO.

# **Additional Considerations**

CT26.WT-iRFP-Neo cells easily detach from tissue culture surfaces. Coating culture plates with poly-D-Lysine prior to use can be used to increase cell adherence to the plates if necessary.

# CT26.WT-iRFP-Neo




# **Certificate of Analysis**

Testing performed by Imanis Life Sciences


Test description	Result
Post thaw viable cell recovery	Pass QC
Sterility	No contamination detected
Mycoplasma	No contamination detected
Neomycin selection	Pass QC
Fluorescence expression	Pass QC

# Morphology:



Low and high density photos taken 24 and 72 hours after thawing, respectively.

# Fluorescence Expression



CT26.WT-iRFP-Neo cells (red) or isotype control (CT26.WT-Fluc-Neo; grey) cells were fixed with paraformaldehyde and analyzed by flow cytometry (20,000 events).

# **Legal Disclaimers**

Limited Product Warranty

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including, without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Imanis. Imanis shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use

For in vitro use only. This certificate is a declaration of analysis at the time of manufacture.

#### PURCHASER NOTIFICATION

#### LIMITED LICENSE NOTICE - RESEARCH USE ONLY

IMANIS LIFE SCIENCES HAS A LIMTED LICENSE UNDER PATENTS OWNED BY THE SALK INSTITUTE FOR BIOLOGICAL STUDIES THAT PERMITS IMANIS LIFE SCIENCES TO SELL PRODUCTS CONTAINING WPRE FOR RESEARCH USE ONLY AND NOT FOR ANY COMMERCIAL USES. EXCLUDED COMMERCIAL USES INCLUDE WITHOUT LIMITATION MANUFACTURING, PROVIDING A SERVICE, THERAPEUTIC, DIAGNOSTIC AND PROPHYLACTIC USES, AND ANY OTHER COMMERCIAL USES. USE OF THIS PRODUCT BY A PURCHASER FOR ANY PURPOSE OTHER THAN FOR RESEARCH IS UNAUTHORIZED AND PROHIBITED.

The Salk Institute actively licenses its patents for commercial uses, and a commercial use license may be available for Salk's WPRE patents. If you wish to inquire about such a license, please

Office of Technology Development The Salk Institute for Biological Studies 10010 North Torrey Pines Road La Jolla, CA 92037 Phone: (858) 453-4100 extension 1278 Fax: (858) 546-8093

Quality control by: RLV Quality Assurance by: SPR Effective Date: 12/07/15