A549-hNIS-Neo



Product Description

Product Name: A549-hNIS-Neo

Catalog Number: CL004 Lot Number: CL-MA04

Species: Human (Homo sapiens)

Strain: Lung

Cell type: Adenocarcinoma

Parental cells: A549 (ATCC® CCL-185TM)*

Morphology: Epithelial Growth mode: Adherent

Reporter gene: Human sodium iodide symporter (hNIS)

Selection gene: Neomycin (Neo)

A549-hNIS-Neo is a polyclonal population of the human lung carcinoma A549 cell line transduced with a lentiviral vector (Imanis #LV013) encoding the human sodium iodide symporter (hNIS) cDNA under the spleen focus-forming virus (SFFV) promoter linked to the neomycin resistance gene via an IRES. High NIS expressing cells were selected using G418. The lentiviral vectors are self-inactivating (SIN) vectors in which the viral enhancer and promoter have been deleted. Transcription inactivation of the LTR in the SIN provirus increases biosafety by preventing mobilization by replication competent viruses and enables regulated expression of the genes from the internal promoters without *cis*-acting effects of the LTR¹.

Mycoplasma Testing

This cell line has tested negative for mycoplasma contamination.

Cell line Authentication

Authentication of the parental A549 cell line was confirmed by short tandem repeat (STR) profiling.

Recommended Uses

In vitro: This is a high hNIS expressing clone suitable for use as a positive control cell line in I-125 uptake assays to validate NIS expression in your lentiviral transduced cells.

In vivo: A549 cells form metastases in the lungs of mice post systemic administration. The in vivo growth of these metastases can be monitored noninvasively in animals using SPECT or PET imaging

References

¹Miyoshi et al. J Virol. 1998. 72:8150-8157.

Storage Instructions

Remove cells from the dry ice packaging and immediately store in the vapor phase above liquid nitrogen (below -130°C).

Complete Growth Medium

Dulbecco's Modified Eagle's Medium (DMEM) 10% Fetal bovine serum (FBS) 1% Penicillin/Streptomycin 0.6 mg/mL G418

G418 should <u>NOT</u> be added to the medium until a culture has been well established from the thawed cells (about 1 week). It is also recommended that a backup frozen cell stock be generated (see below) before adding G418 to the growth medium.

Thawing Instructions

- 1. Thaw cells by gently swirling in a 37°C water bath. To limit contamination, do not submerge the O-ring and cap.
- 2. When cells are ~70% thawed (~1 min), remove the vial and wipe down with 70% ethanol. Allow tube to dry completely.
- In a biosafety cabinet, transfer the cells into a 15 mL conical tube containing 5 mL of pre-warmed complete growth medium. Centrifuge cells at ~250 x g for 3-5 min.
- Remove supernatant and resuspend cells in 1 mL complete growth medium. Transfer cells to a T75 flask containing 10 mL pre-warmed complete growth medium.
- 5. Incubate the culture at 37°C with 5% CO₂. Cells should reach full confluency 1-2 days after thawing.

Subculturing Instructions

Volumes are given for a T75 flask. Increase or decrease as needed.

- 1. Remove culture medium from cells.
- Carefully wash the cell monolayer with 5-10 mL of phosphate buffered saline.
- Add 2 mL of 0.25% Trypsin-EDTA solution to the flask and incubate at room temperature until cells have dissociated (approx. 2-5 min).
- 4. Neutralize the trypsin by adding 8 mL complete growth medium, and mix by gently pipetting up and down.
- 5. Transfer desired portion of the cells to a fresh T75 flask. Add fresh complete growth medium to a total volume of 10 mL and return cells to 37°C/5% CO₂ incubator.

For maintenance, a subcultivation ratio of 1:10 is recommended. At this ratio cells will be ready for passage approximately every 3-4 days.

Freezing Medium

Cells can be amplified and used to generate additional frozen stocks. Frozen stocks should be preserved in a designated cryopreservation medium or in complete growth medium <u>without G418</u> supplemented with 5-10% DMSO.

^{*} The ATCC trademark and any and all ATCC catalog numbers are trademarks of the American Type Culture Collection

A549-hNIS-Neo





Certificate of Analysis

Testing performed by Imanis Life Sciences:

Test description	Result
Post thaw viable cell recovery	95%
Sterility	No contamination detected
Mycoplasma	No contamination detected
Neomycin selection	Pass QC
125I uptake assay	Pass QC


Morphology:





Low and high density photos taken at various times.

125 luptake:



Uptake of 125 by 2 x 105 cells was assayed in the presence or absence of KCIO₄, an inhibitor of NIS-mediated ¹²⁵I uptake.

Legal Disclaimers

LIMITED PRODUCT WARRANTY

THIS WARRANTY LIMITS OUR LIABILITY TO REPLACEMENT OF THIS PRODUCT. NO OTHER WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE PROVIDED BY IMANIS. IMANIS SHALL HAVE NO LIABILITY FOR ANY DIRECT, INDIRECT, CONSEQUENTIAL, OR INCIDENTAL DAMAGES ARISING OUT OF THE USE, THE RESULTS OF USE, OR THE INABILITY TO USE THIS PRODUCT.

FOR IN VITRO USE ONLY. THIS CERTIFICATE IS A DECLARATION OF ANALYSIS AT THE TIME OF MANUFACTURE.

LIMITED LICENSE NOTICE - RESEARCH USE ONLY

IMANIS LIFE SCIENCES HAS A LIMTED LICENSE UNDER PATENTS OWNED BY THE SALK INSTITUTE FOR BIOLOGICAL STUDIES THAT PERMITS IMANIS LIFE SCIENCES TO SELL PRODUCTS CONTAINING WPRE FOR RESEARCH USE ONLY AND NOT FOR ANY COMMERCIAL USES. EXCLUDE COMMERCIAL USES INCLUDE WITHOUT LIMITATION MANUFACTURING, PROVIDING A SERVICE, THERAPEUTIC, DIAGNOSTIC AND PROPHYLACTIC USES, AND ANY OTHER COMMERCIAL USES. USE OF THIS PRODUCT BY A PURCHASER FOR ANY PURPOSE OTHER THAN FOR RESEARCH IS

THE IMANIS CELL LINES ARE NOT INTENDED FOR USE IN HUMANS. CELL LINES TRANSDUCED WITH LENTIVIRAL VECTORS ARE CLASSIFIED AS BIOSAFETY LEVEL 2 REAGENTS AND SHOULD BE USED UNDER THE APPROPRIATE

THE IMANIS CELL LINES ARE NOT INTENDED FOR USE IN HUMANS. CELL LINES TRANSDUCED WITH LENTIVIRAL VECTORS ARE CLASSIFIED AS BIOSAFETY LEVEL PER REAGENTS AND SHOULD BE USED UNDER THE APPROPRIATE BIOSAFETY LEVEL PER INSTITUTIONAL GUIDELINES.

THE PURCHASER AGREES THAT IMANIS MATERIALS DESIGNATED AS BIO-SAFETY LEVEL 2 CONSTITUTE KNOWN PATHOGENS AND THAT OTHER IMANIS MATERIALS NOT SO DESIGNATED AND ANY PROGENY OR MODIFICATION MAY BE PATHOGENIC UNDER CERTAIN CONDITIONS. PURCHASER ASSUMES ALL RISK AND RESPONSIBILITY IN CONNECTION WITH THE RECEIPT, HANDLING, STORAGE, DISPOSAL, TRANSFER AND USE OF THE IMANIS MATERIALS INCLUDING WITHOUT LIMITATION TAKING ALL APPROPRIATE SAFETY AND HANDLING PRECAUTIONS TO MINIMIZE HEALTH OR ENVIRONMENTAL RISK, PURCHASER AGREES THAT ANY ACTIVITY UNDERTRAKEN WITH THE IMANIS MATERIALS AND ANY PROGENY OR MODIFICATION WILL BE CONDUCTED IN COMPLIANCE WITH ALL APPLICABLE GUIDELINES, LAWS AND REGULATIONS. THE IMANIS MATERIALS ARY OTHER IMANIS PROJECT ON THE MANIS MATERIALS AND ANY PROGENY OR MODIFICATION WITH DET OR ANY MITHOUT WARRANTIES OF MAY TECHNICAL INFORMATION AND ASSISTANCE PROVIDED BY IMANIS ARE PROVIDED AS IS, WITHOUT WARRANTIES OF MY KIND, EXPRESS OR IMPLIED, INCLUDINGS BUT NOT LIMITED TO ANY IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, MANUFACTURE ACCORDING TO GMP STANDARDS, TYPICALITY, SAFETY, ACCURACY AND NON-INFRINGEMENT. IN NO EVENT SHALL IMANIS, ITS PARENTS, SUBSIDIARIES, DIRECTORS, OFFICERS, AGENTS, EMPLOYES, ASSIGNS, SUCCESSORS AND AFFILIATE (COLLECTIVELY "IMANIS INDEMNIFIED PARTIES) BE LIBLE FOR INDIRECT, SPECIAL, INCIDENTAL OR CONSCIUNTAL LAWAS. ITS PARENTS, SUBSIDIARIES DIRECTORS, OFFICERS, AGENTS, EMPLOYES, ASSIGNS, SUCCESSORS AND AFFILIATE (COLLECTIVELY "IMANIS INDEMNIFIED PARTIES) BE LIBLE FOR INDIRECT, SPECIAL, INCIDENTAL OR CONSCIUNTAL IMANIS OF THE POSSIBILITY OF SUCH DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOST PROFITS, COST OF CAPITAL, COST OF SUBSTITUTE PRODUCTS OR CLAIMS OF LICENSEES CUSTOMERS FOR SUCH DAMAGES. IN CUEVET SHALL IMANIS CITI

Quality control by: RLV/AWD Quality Assurance by: RLV/SPR Effective Date: 20-Mar-2024