## **Vesicular Stomatitis Virus (VSV) antibody**



#### **Product Description**

Product Name: Vesicular Stomatitis Virus (VSV) antibody

Catalog Number: REA005
Lot Number: IM04
Unit Size: 250 μL
Species: Rabbit
Antigen: VSV-hIFNβ
Isotype: Polyclonal
Reactivity: VSV

This is a rabbit polyclonal antibody is raised against a recombinant Indiana strain of Vesicular Stomatitis Virus (VSV) encoding human interferon beta (VSV-hIFN $\beta$ ). The antibody detects the VSV G, N, and M proteins.

#### **Storage Instructions**

This antibody can be stored short term (1-2 weeks) at 2-6°C. For longer term storage, aliquot and store at or below -20°C. Avoid repeated freeze/thaw cycles.

#### **Recommended Uses**

The VSV antibody is suitable for a number of different applications, including:

Immunoblot (1:8000-1:10,000 dilution)
Immunofluorescence (1:1500-1:3000 dilution)
Immunohistochemistry (1:5000-1:10,000 dilution)

Flow cytometry (1:10,000 dilution)

ELISA (1:2000-1:4000 dilution to detect 1 x  $10^7$  TCID<sub>50</sub>)


Antibody neutralization (see below)

The VSV antibody has been tested for neutralization against VSV-GFP. The neutralizing titer is 1:2560 for complete inhibition of 400 TCID<sub>50</sub> units of VSV-GFP on Vero cells.

#### **Certificate of Analysis**


Testing performed by Imanis Life Sciences

#### **Immunoblot**



Lysates prepared from VSV-GFP-infected (lane 1) Vero cells or VSV-GFP particles (lane 2) were subjected to SDS-PAGE/ immunoblot using anti-VSV antibody (1:10,000 dilution) and HRP-conjugated anti-rabbit secondary antibody.


#### **Immunofluorescence**





VSV-mIFNβ-NIS- (left) or saline- (right) treated tumor xenografts were stained with VSV antibody (1:1500 dilution) followed by an Alexa Fluor 594-conjugated anti rabbit secondary antibody.


#### **Immunohistochemistry**





VSV-mIFNβ-NIS- (left) or saline- (right) treated tumor xenografts were stained with VSV antibody (1:5000 dilution). Counterstain is hematoxylin.

## Flow Cytometry



VSV-GFP-infected Vero cells were fixed with paraformaldehyde and stained with prebleed antiserum (red) or anti-VSV antibody (blue; dilution 1:10,000) followed by an Alexa Fluor 594-conjugated secondary antibody. The stained cells were subjected to flow cytometry analysis.

# **Vesicular Stomatitis Virus (VSV) antibody**



### **Legal Disclaimers**

LIMITED PRODUCT WARRANTY
THIS WARRANTY LIMITS OUR LIABILITY TO REPLACEMENT OF THIS PRODUCT. NO OTHER WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE PROVIDED BY IMANIS. IMANIS SHALL HAVE NO LIABILITY FOR ANY DIRECT, INDIRECT, INDIR CONSEQUENTIAL, OR INCIDENTAL DAMAGES ARISING OUT OF THE USE, THE RESULTS OF USE, OR THE INABILITY TO USE THIS PRODUCT.

FOR IN VITRO USE ONLY. THIS CERTIFICATE IS A DECLARATION OF ANALYSIS AT THE TIME OF MANUFACTURE.

#### PURCHASER NOTIFICATION

#### LIMITED LICENSE NOTICE - RESEARCH USE ONLY

THE NIS GENE AND TECHNOLOGY IS COVERED UNDER AN EXCLUSIVE LICENSE TO IMANIS LIFE SCIENCES. RESEARCHERS MAY USE THIS PRODUCT FOR RESEARCH USE ONLY AND LIFE SCIENCES. RESEARCHERS MAY USE I HIS PRODUCT FOR RESEARCH USE ONLY AND NO COMMERCIAL USE IS ALLOWED. NO OTHER USE OR TRANSFER OF THIS PRODUCT OR DERIVATIVES IS AUTHORIZED WITHOUT THE PRIOR EXPRESS WRITTEN CONSENT OF IMANIS. WITH RESPECT TO ANY USES OUTSIDE THIS LABEL LICENSE, INCLUDING ANY DIAGNOSTIC, THERAPEUTIC OR PROPHYLACTIC USES, PLEASE CONTACT IMANIS LIFE SCIENCES FOR SUPPLY AND LICENSING INFORMATION.

Imanis Life Sciences 2900 37<sup>th</sup> Street NW, Building 110 Rochester MN 55901 Tel: (507) 218-2559 Email: support@imanislife.com